



A New Addition  
to the **MEDRAY**  
Family of Lasers

**SUPERMAX**  
MEDRAY CLASS IV LASER

MEDRAY SUPERMAX is one of the strongest lasers on the market and is combined with an optimal set of innovative features.

## FOUR WAVELENGTHS: 650 | 810 | 915 | 1064

These waves are the absolute best for, not only **deep penetration and healing**, but the utmost in **pain control**, as well as **photo-biomodulation**.

MEDRAY SUPERMAX upholds the MEDRAY Class IV Laser practitioner promise of delivering the best in:

 **QUALITY**     **SERVICE**     **AFFORDABILITY**

### ADDITIONAL KEY SPECIFICATIONS

**Laser Type:** GaAlAs Diode Laser System | Class IV

**Max Power:** 300mw+15W+15W+15W (Total 45.3W)

**Working Mode:** Pulsed or Continuous Wave (CW)

**Dimensions & Weight:** 275(L)\*283(W)\*186(H) mm | 7Kg



# A PERFECT BLEND OF POWER, INNOVATION, & AFFORDABILITY

## CUSTOM WAVELENGTH PROTOCOLS

Each wavelength (650 | 810 | 915 | 1064) can deliver 15 watts, making it powerful enough to handle extremely challenging and/or complex cases.

Create and save custom protocols for specific patients, making treatment more efficient than ever before! Additionally, there are 16 preset protocols based for the areas of pain most commonly treated.

## OUTSTANDING ACCESSORIES

Providing patients with the most effective and efficient laser treatment has never been easier with the accessories offered with the MEDRAY SUPERMAX.

You can use the standard, attached emitter head, or quickly remove it for the brand new, state-of-the-art, roller headpiece. You can also be totally hands-free with the patented, optional S.T.A.R. System. S.T.A.R. includes an elevated cart, custom flex arm, 3D emitter holder, wireless stop switch and optical headpiece.



## INCREDIBLE GRAPHICS

The graphics provided on the MEDRAY SUPERMAX exude the “wow” factor patients seek from modern practices.

**STAR.**  
STATIONARY THERAPEUTIC APPLICATION OF RADIANCE

## AFFORDABLE PRICING

We are known for offering an unmatched combination of technology and pricing.

MEDRAY SUPERMAX: \$16,995  
OPTIONAL S.T.A.R. \$1,495

**SUPERMAX**  
MEDRAY CLASS IV LASER



Contact Dr. Gary Huddleston to learn more about providing your patients with the elite technology of MEDRAY SUPERMAX | 573.745.1086 | DrG4000@att.net | [www.medraylaser.com](http://www.medraylaser.com)